

# Anti-Myc Magnetic Beads

## Product Description

The Biolinkedin® Anti-Myc Magnetic Beads are used for immunoprecipitation (IP) of specific Myc-tagged proteins expressed in human in vitro expression systems and bacterial and mammalian cell lysates. Anti-Myc magnetic beads are based on carboxyl magnetic beads, with 200 nm particle size, covalently coupling with high quality mouse IgG monoclonal antibody that recognizes the Myc-epitope tag (EQKLISEEDL) derived from the human Myc protein. For IP, the beads are added to a sample containing Myc-tagged proteins. The bound Myc-tagged proteins are dissociated from the beads using an elution buffer. The beads are removed from the solution manually using a magnetic stand or by automation

## Product Features

Composition	mouse IgG <sub>1</sub> monoclonal Ab
Magnetization	Superparamagnetic
Particle size	200 nm
Concentration	10 mg/mL
Binding Capacity	≥ 0.5 mg Myc-tagged fusion protein/mL of beads
Application	IP, CoIP
Storage Condition	Store at 4°C for 2 years.

## Protocol

### 1. Cell lysis

Cells may be lysed using any standard cell lysis protocol compatible with your starting material. We recommend the use of **Biolinkedin® IP Lysis/Wash Buffer**. Add protease inhibitor (such as PMSF at 1mM) if needed.

### 2. Preparation of Magnetic Beads

2.1 Resuspend the Magnetic Beads in the vial (tilt and rotate for 2 minutes or gently pipette for 10 times).

2.2 Transfer 25-50 µL of Anti-Myc Magnetic Beads into a 1.5 mL tube (Transfer amount may be adjusted as required).

2.3 Add 500 µL of IP Lysis/Wash Buffer to the beads and gently pipette to mix. Place the tube into a magnetic stand to collect the beads against the side of the tube (Herein after referred to as magnetic separation).

Remove and discard the supernatant. Repeat this step twice.

### 3. Immunoprecipitation

3.1 Remove the tubes from the magnetic separator and add your sample containing Myc-tagged protein to the pre-washed magnetic beads and incubate at room temperature for 2h or overnight at 4°C with mixing.

3.2 Collect the beads with a magnetic stand, remove the unbound sample and collect for analysis.

3.3 Add 500µL of IP Lysis/Wash Buffer to the tube and gently mix. Collect the beads and discard the supernatant. Repeat wash twice.

### 4. Elution

**Note:** Select one of the elution protocols below. If the eluted Myc-tagged protein will be used for function applications or is sensitive to pH extremes, then elute the protein with HA Peptide.

#### Gentle Elution Protocol

1. Prepare Myc Peptide at 0.2-1mg/mL in TBS.
2. Add 100 µL of 0.2-1mg/mL Myc Peptide to the beads, gently vortex to mix and incubate the sample at 4°C on a rotator for 2h-6h.
3. Separate the beads on a magnetic stand and save the supernatant containing the target antigen.
4. Repeat elution step once for higher recovery.

#### Chemical Elution Protocols

1. Add 100µL of 1x SDS-PAGE loading buffer to the tube.
2. Boil for 5 minutes on a dry bath.
3. Magnetically separate the beads and save the supernatant containing the target antigen.
  - Acidic Elution
    1. Add 100µL of 0.1M glycine, pH 3.0.
    2. Gently vortex to mix and incubate the sample at room temperature on a rotator for 5-10 minutes.
    3. Magnetically separate the beads and save the supernatant containing the target antigen.
    4. To neutralize the low pH, add 20µL of Neutralization Buffer (1 M Tris pH 8.5) for each 100µL of eluate.

The final solution can be used as samples for denaturing SDS-PAGE. Or the elution can be adjusted to neutral pH with neutralization buffer immediately and used for further analysis.

## Troubleshooting

Problem	Possible Cause	Solution
High background band	Non-specific binding of proteins to antibodies, insufficient washing on magnetic beads or EP tubes	The lysate was pretreated to remove non-specific proteins; before the last wash, the entire sample was transferred to a new EP tube and centrifuged.
	Not enough washing	Increase the time and number of washes
No protein band	No or minimal tagged protein was expressed	Verify protein expression; Prepare fresh lysate; Use the appropriate protease inhibitor.
	Insufficient incubation time	Increase incubation time
	Interfering substances in the sample	High concentrations of DTT, 2-mercaptoethanol or other reducing agents are present in the lysate.