Biolinkedin® 基础 His 标签蛋白纯化试剂盒(缓冲液)

包装清单

名称	货号	包装
基础 His 标签蛋白	BK-2009	<u>~</u>
纯化试剂盒 (缓冲液)	BK-2009	盒

产品概述

Biolinkedin[®] 基础 His 标签蛋白纯化试剂盒(缓冲液) 是由优化预制的缓冲液组成,可搭配 His 标签蛋白琼脂糖 凝胶或 His 标签蛋白琼脂糖磁珠,用于纯化各种表达系统 融合表达的 His 标签重组蛋白。

试剂盒组成

结合/平衡缓冲液 (2×)	100 mL
洗涤缓冲液 (5×)	100 mL
洗脱缓冲液	100 mL

操作流程

操作流程可参考 PK-2001 His 标签蛋白纯化试剂盒(凝胶法)、PK-2003 His 标签蛋白纯化试剂盒或 PK-2007 His 标签蛋白纯化试剂盒(NTA-Ni 凝胶法)说明书具体操作步骤。

相关产品

THIN HH	
货号	产品名称
PK-2001	His 标签蛋白纯化试剂盒(凝胶法)
PK-2003	His 标签蛋白纯化试剂盒
PK-2005	His 标签蛋白纯化再生试剂盒
PK-2007	His 标签蛋白纯化试剂盒 (NTA-Ni 凝胶法)
BK-2008	His 标签蛋白纯化试剂盒(包涵体蛋白)
BK-2009	基础 His 标签蛋白纯化试剂盒 (缓冲液)